
1 SIMPATI* TCP/IP-Interface (simserv)

1.1 Description

Other applications like LabView, HP VEE Lab or Siemens-S7 can communicate via **the TCP/IP socket interface** from Windows with SIMPATI and can exchange all relevant process data.


Hint:

The TCP/IP-interface will be supported from SIMPATI version 2.0 on. In older versions maybe not all commands available.


1.2 Function

Attached the function is pictured:

1.3 Installation

The TCP/IP-Server software will be started with the entry **80:01:simserv:-P7777:** in the file **simpati.str** in the directory **..\simpati\system**. After that a new start of **SIMPATI** is necessary.

Close **SIMPATI** before you change the configuration. Open the file with the windows standard-program Notepad/Editor and use a continuous numbering for the entries! Incorrect entries can lead to a crash of the program **SIMPATI**.

1.3.1 Command procedure

It is not possible to open a permanent connection to the server. Therefore a server can operate with a lot of clients.

For each command the following terminate is to be used:

1. Open connection
2. Send command
3. Read response
4. Close connection

1.4 Function-commands

Attached find the description of the function-commands (Simserv-commands for the data exchange with [SIMPATI](#)).

1.4.1 Command Syntax

The command syntax is as follows:

| Syntax | | | | | | | | | |
|---------|----|------------|----|----------|----|----------|----|--------------------------|----|
| Command | SR | Chamber-ID | SR | Param. 1 | SR | Param. 2 | SR | Parameter n (up to 4) | CR |

SR = Separator = "¶" (ASCII 182)

CR = CarriageReturn (\r) (ASCII 13)

In the event of an error:

| | |
|-------------------|-----------|
| Error code | CR |
|-------------------|-----------|

Possible error codes:

| Error code | Error text |
|------------|---|
| -1 | Empty command string |
| -2 | Missing Chamber-ID |
| -3 | Invalid Chamber-ID |
| -4 | Chamber not accessible |
| -5 | Unknown Command-ID |
| -6 | Insufficient number or wrong command parameters |
| -7 | No server (in server service command mode) |

In a fault-free case a "1" is returned.

1.4.2 Examples

Set setpoint Temperature of the 2nd Chamber to 25°C:

| | | | | | | | |
|-------|---|---|---|---|---|------|----|
| 11001 | ¶ | 2 | ¶ | 1 | ¶ | 25.0 | CR |
|-------|---|---|---|---|---|------|----|

Hint:

For the first control value Temperature usually the value 1 and for the second control value (e.g. Humidity) the value 2 is used etc.!

Set actual value Temperature of the 2nd Chamber to 30°C:

| | | | | | | | |
|-------|---|---|---|---|---|------|----|
| 11003 | ¶ | 2 | ¶ | 1 | ¶ | 30.0 | CR |
|-------|---|---|---|---|---|------|----|

Read the actual value Temperature of the 2nd Chamber:

| | | | | | | |
|-------|---|---|---|---|---|----|
| 11004 | ¶ | 2 | ¶ | 1 | ¶ | CR |
|-------|---|---|---|---|---|----|

Response from the server on a command GET ISTWERT:

| | | | |
|---|---|---------|----|
| 1 | ¶ | 23.9000 | CR |
|---|---|---------|----|

Read the digital channel 1(Start) of the 2nd Chamber:


| | | | | | | |
|-------|---|---|---|---|---|----|
| 14003 | ¶ | 2 | ¶ | 1 | ¶ | CR |
|-------|---|---|---|---|---|----|

Set the digital channel 1(Start) of the 2nd Chamber:

| | | | | | | | |
|-------|---|---|---|---|---|---|----|
| 14001 | ¶ | 2 | ¶ | 1 | ¶ | 1 | CR |
|-------|---|---|---|---|---|---|----|


1.5 Test program

The TCP/IP connection can be checked with the test program „**Servtest.exe**“, which is located in the directory `..\simpati\system`.
(before execute point 1.3)

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